



NUP93 Antibody, FITC conjugated

Product Code	CSB-PA822683LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N1F7
Immunogen	Recombinant Human Nuclear pore complex protein Nup93 protein (14-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup93 (93 kDa nucleoporin) (Nucleoporin Nup93), NUP93, KIAA0095
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NUP93
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.