



EPS8L3 Antibody, HRP conjugated

Product Code	CSB-PA823469LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE67
Immunogen	Recombinant Human Epidermal growth factor receptor kinase substrate 8-like protein 3 protein (367-521AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epidermal growth factor receptor kinase substrate 8-like protein 3 (EPS8-like protein 3) (Epidermal growth factor receptor pathway substrate 8-related protein 3) (EPS8-related protein 3), EPS8L3, EPS8R3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EPS8L3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.