



EPS8L3 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA823469LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8TE67 |
| Immunogen | Recombinant Human Epidermal growth factor receptor kinase substrate 8-like protein 3 protein (367-521AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Epidermal growth factor receptor kinase substrate 8-like protein 3 (EPS8-like protein 3) (Epidermal growth factor receptor pathway substrate 8-related protein 3) (EPS8-related protein 3), EPS8L3, EPS8R3 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | EPS8L3 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |