



SLC45A3 Antibody, HRP conjugated

Product Code	CSB-PA836239LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JT2
Immunogen	Recombinant Human Solute carrier family 45 member 3 protein (408-487AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Solute carrier family 45 member 3 (Prostate cancer-associated protein 6) (Prostein), SLC45A3, PCANAP6 PRST
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SLC45A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.