



POFUT2 Antibody, FITC conjugated

Product Code	CSB-PA837737LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2G5
Immunogen	Recombinant Human GDP-fucose protein O-fucosyltransferase 2 protein (125-271AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GDP-fucose protein O-fucosyltransferase 2 (EC 2.4.1.221) (Peptide-O-fucosyltransferase 2) (O-FucT-2), POFUT2, C21orf80 FUT13 KIAA0958
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	POFUT2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.