



CNNM3 Antibody, HRP conjugated

Product Code	CSB-PA843333LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE01
Immunogen	Recombinant Human Metal transporter CNNM3 protein (358-645AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metal transporter CNNM3 (Ancient conserved domain-containing protein 3) (Cyclin-M3), CNNM3, ACDP3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CNNM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.