



MASTL Antibody, HRP conjugated

Product Code	CSB-PA846629LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96GX5
Immunogen	Recombinant Human Serine/threonine-protein kinase greatwall protein (612-740AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase greatwall (GW) (GWL) (hGWL) (EC 2.7.11.1) (Microtubule-associated serine/threonine-protein kinase-like) (MAST-L), MASTL, GW GWL THC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MASTL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.