



ATL1 Antibody, Biotin conjugated

Product Code	CSB-PA855500LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXF7
Immunogen	Recombinant Human Atlastin-1 protein (220-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Atlastin-1 (EC 3.6.5.-) (Brain-specific GTP-binding protein) (GTP-binding protein 3) (GBP-3) (hGBP3) (Guanine nucleotide-binding protein 3) (Spastic paraplegia 3 protein A), ATL1, GBP3 SPG3A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ATL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.