



DVL3 Antibody, HRP conjugated

Product Code	CSB-PA856442LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92997
Immunogen	Recombinant Human Segment polarity protein dishevelled homolog DVL-3 protein (509-696AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Segment polarity protein dishevelled homolog DVL-3 (Dishevelled-3) (DSH homolog 3), DVL3, KIAA0208
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DVL3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.