



LPP Antibody, HRP conjugated

Product Code	CSB-PA856611LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q93052
Immunogen	Recombinant Human Lipoma-preferred partner protein (227-391AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipoma-preferred partner (LIM domain-containing preferred translocation partner in lipoma), LPP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LPP