



DUSP9 Antibody, Biotin conjugated

Product Code	CSB-PA859528LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99956
Immunogen	Recombinant Human Dual specificity protein phosphatase 9 protein (142-263AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity protein phosphatase 9 (EC 3.1.3.16) (EC 3.1.3.48) (Mitogen-activated protein kinase phosphatase 4) (MAP kinase phosphatase 4) (MKP-4), DUSP9, MKP4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DUSP9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.