



MYOC Antibody, HRP conjugated

Product Code	CSB-PA859950LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99972
Immunogen	Recombinant Human Myocilin protein (183-294AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myocilin (Myocilin 55 kDa subunit) (Trabecular meshwork-induced glucocorticoid response protein) [Cleaved into: Myocilin, N-terminal fragment (Myocilin 20 kDa N-terminal fragment); Myocilin, C-terminal fragment (Myocilin 35 kDa N-terminal fragment)], MYOC, GLC1A TIGR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MYOC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.