



ABHD6 Antibody, Biotin conjugated

Product Code	CSB-PA861144LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV23
Immunogen	Recombinant Human Monoacylglycerol lipase ABHD6 protein (105-211AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Monoacylglycerol lipase ABHD6 (EC 3.1.1.23) (2-arachidonoylglycerol hydrolase) (Abhydrolase domain-containing protein 6), ABHD6
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ABHD6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.