



CLEC4M Antibody, Biotin conjugated

Product Code	CSB-PA862003LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H2X3
Immunogen	Recombinant Human C-type lectin domain family 4 member M protein (238-382AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 4 member M (CD209 antigen-like protein 1) (DC-SIGN-related protein) (DC-SIGNR) (Dendritic cell-specific ICAM-3-grabbing non-integrin 2) (DC-SIGN2) (Liver/lymph node-specific ICAM-3-grabbing non-integrin) (L-SIGN) (CD antigen CD299), CLEC4M, CD209L CD209L1 CD299
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CLEC4M
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.