



NAT10 Antibody, FITC conjugated

Product Code	CSB-PA863939LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H0A0
Immunogen	Recombinant Human RNA cytidine acetyltransferase protein (909-1025AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA cytidine acetyltransferase (EC 2.3.1.-) (18S rRNA cytosine acetyltransferase) (N-acetyltransferase 10), NAT10, ALP KIAA1709
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAT10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.