



NR2E3 Antibody, HRP conjugated

Product Code	CSB-PA865585LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5X4
Immunogen	Recombinant Human Photoreceptor-specific nuclear receptor protein (112-222AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Photoreceptor-specific nuclear receptor (Nuclear receptor subfamily 2 group E member 3) (Retina-specific nuclear receptor), NR2E3, PNR RNR
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NR2E3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.