



ABCD2 Antibody, Biotin conjugated

Product Code	CSB-PA866202LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UBJ2
Immunogen	Recombinant Human ATP-binding cassette sub-family D member 2 protein (559-715AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family D member 2 (Adrenoleukodystrophy-like 1) (Adrenoleukodystrophy-related protein) (hALDR), ABCD2, ALD1 ALDL1 ALDR ALDRP
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ABCD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.