



RNF19A Antibody, Biotin conjugated

Product Code	CSB-PA868327LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NV58
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF19A protein (741-825AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF19A (EC 2.3.2.27) (Double ring-finger protein) (Dorfin) (RING finger protein 19A) (p38), RNF19A, RNF19
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RNF19A