



MID1IP1 Antibody, Biotin conjugated

Product Code	CSB-PA873607LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPA3
Immunogen	Recombinant Human Mid1-interacting protein 1 protein (1-183AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mid1-interacting protein 1 (Gastrulation-specific G12-like protein) (Mid1-interacting G12-like protein) (Protein STRAIT11499) (Spot 14-related protein) (S14R) (Spot 14-R), MID1IP1, MIG12
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MID1IP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.