



RHCG Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA880059LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UBD6 |
| Immunogen | Recombinant Human Ammonium transporter Rh type C protein (416-479AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Ammonium transporter Rh type C (Rh glycoprotein kidney) (Rhesus blood group family type C glycoprotein) (Rh family type C glycoprotein) (Rh type C glycoprotein) (Tumor-related protein DRC2), RHCG, C15orf6 CDRC2 PDRC2 RHGK |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | RHCG |