



MMP28 Antibody, FITC conjugated

Product Code	CSB-PA880950LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H239
Immunogen	Recombinant Human Matrix metalloproteinase-28 protein (351-475AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Matrix metalloproteinase-28 (MMP-28) (EC 3.4.24.-) (Epilysin), MMP28, MMP25
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MMP28
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.