



FASTKD2 Antibody, FITC conjugated

Product Code	CSB-PA882149LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYY8
Immunogen	Recombinant Human FAST kinase domain-containing protein 2, mitochondrial protein (117-203AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FAST kinase domain-containing protein 2, mitochondrial, FASTKD2, KIAA0971
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	FASTKD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.