



MKRN1 Antibody, Biotin conjugated

Product Code	CSB-PA883390LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UHC7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase makorin-1 protein (109-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase makorin-1 (EC 2.3.2.27) (RING finger protein 61) (RING-type E3 ubiquitin transferase makorin-1), MKRN1, RNF61
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MKRN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.