



LGR4 Antibody, HRP conjugated

Product Code	CSB-PA883618LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXB1
Immunogen	Recombinant Human Leucine-rich repeat-containing G-protein coupled receptor 4 protein (234-365AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat-containing G-protein coupled receptor 4 (G-protein coupled receptor 48), LGR4, GPR48
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LGR4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.