



ULBP1 Antibody, HRP conjugated

Product Code	CSB-PA887177LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZM6
Immunogen	Recombinant Human UL16-binding protein 1 protein (26-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UL16-binding protein 1 (ALCAN-beta) (NKG2D ligand 1) (N2DL-1) (NKG2DL1) (Retinoic acid early transcript 1I), ULBP1, N2DL1 RAET1I
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ULBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.