



# ULBP1 Antibody, Biotin conjugated

Product Code	CSB-PA887177LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZM6
Immunogen	Recombinant Human UL16-binding protein 1 protein (26-216AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UL16-binding protein 1 (ALCAN-beta) (NKG2D ligand 1) (N2DL-1) (NKG2DL1) (Retinoic acid early transcript 1I), ULBP1, N2DL1 RAET1I
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ULBP1