



HRH4 Antibody, HRP conjugated

Product Code	CSB-PA887972LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3N8
Immunogen	Recombinant Human Histamine H4 receptor protein (204-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Histamine H4 receptor (H4R) (HH4R) (AXOR35) (G-protein coupled receptor 105) (GPRv53) (Pfi-013) (SP9144), HRH4, GPCR105
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HRH4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.