



WASF3 Antibody, FITC conjugated

Product Code	CSB-PA891467LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPY6
Immunogen	Recombinant Human Wiskott-Aldrich syndrome protein family member 3 protein (368-502AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wiskott-Aldrich syndrome protein family member 3 (WASP family protein member 3) (Protein WAVE-3) (Verprolin homology domain-containing protein 3), WASF3, KIAA0900 SCAR3 WAVE3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASF3