



WASF3 Antibody, Biotin conjugated

Product Code	CSB-PA891467LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPY6
Immunogen	Recombinant Human Wiskott-Aldrich syndrome protein family member 3 protein (368-502AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wiskott-Aldrich syndrome protein family member 3 (WASP family protein member 3) (Protein WAVE-3) (Verprolin homology domain-containing protein 3), WASF3, KIAA0900 SCAR3 WAVE3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASF3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.