



RPS6KA6 Antibody, Biotin conjugated

Product Code	CSB-PA891943LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK32
Immunogen	Recombinant Human Ribosomal protein S6 kinase alpha-6 protein (340-467AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ribosomal protein S6 kinase alpha-6 (S6K-alpha-6) (EC 2.7.11.1) (90 kDa ribosomal protein S6 kinase 6) (p90-RSK 6) (p90RSK6) (Ribosomal S6 kinase 4) (RSK-4) (pp90RSK4), RPS6KA6, RSK4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RPS6KA6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.