



TPX2 Antibody, FITC conjugated

Product Code	CSB-PA892151LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW0
Immunogen	Recombinant Human Targeting protein for Xklp2 protein (160-320AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Targeting protein for Xklp2 (Differentially expressed in cancerous and non-cancerous lung cells 2) (DIL-2) (Hepatocellular carcinoma-associated antigen 519) (Hepatocellular carcinoma-associated antigen 90) (Protein fls353) (Restricted expression proliferation-associated protein 100) (p100), TPX2, C20orf1 C20orf2 DIL2 HCA519
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TPX2