



TPX2 Antibody, Biotin conjugated

Product Code	CSB-PA892151LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULW0
Immunogen	Recombinant Human Targeting protein for Xklp2 protein (160-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Targeting protein for Xklp2 (Differentially expressed in cancerous and non-cancerous lung cells 2) (DIL-2) (Hepatocellular carcinoma-associated antigen 519) (Hepatocellular carcinoma-associated antigen 90) (Protein fls353) (Restricted expression proliferation-associated protein 100) (p100), TPX2, C20orf1 C20orf2 DIL2 HCA519
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TPX2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.