



EDAR Antibody, HRP conjugated

Product Code	CSB-PA892357LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UNE0
Immunogen	Recombinant Human Tumor necrosis factor receptor superfamily member EDAR protein (27-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tumor necrosis factor receptor superfamily member EDAR (Anhidrotic ectodysplasin receptor 1) (Downless homolog) (EDA-A1 receptor) (Ectodermal dysplasia receptor) (Ectodysplasin-A receptor), EDAR, DL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EDAR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.