



# MAN2B2 Antibody, FITC conjugated

|                     |   |
|---------------------|---|
| Product Code        | CSB-PA896489LC01HU  |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| Uniprot No.         | Q9Y2E5  |
| Immunogen           | Recombinant Human Epididymis-specific alpha-mannosidase protein (347-523AA)                                 |
| Raised In           | Rabbit  |
| Species Reactivity  | Human   |
| Form                | Liquid  |
| Conjugate           | FITC  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                            |
| Purification Method | >95%, Protein G purified  |
| Isotype             | IgG   |
| Clonality           | Polyclonal  |
| Alias               | Epididymis-specific alpha-mannosidase (EC 3.2.1.24) (Mannosidase alpha class 2B member 2), MAN2B2, KIAA0935 |
| Immunogen Species   | Homo sapiens (Human)  |
| Research Area       | Others  |
| Target Names        | MAN2B2  |