



MAN2B2 Antibody, Biotin conjugated

Product Code	CSB-PA896489LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2E5
Immunogen	Recombinant Human Epididymis-specific alpha-mannosidase protein (347-523AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epididymis-specific alpha-mannosidase (EC 3.2.1.24) (Mannosidase alpha class 2B member 2), MAN2B2, KIAA0935
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MAN2B2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.