



PTCH2 Antibody, Biotin conjugated

Product Code	CSB-PA896928LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6C5
Immunogen	Recombinant Human Protein patched homolog 2 protein (793-951AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein patched homolog 2 (PTC2), PTCH2
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	PTCH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.