



SSX2IP Antibody, Biotin conjugated

Product Code	CSB-PA897526LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2D8
Immunogen	Recombinant Human Afadin- and alpha-actinin-binding protein (485-614AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Afadin- and alpha-actinin-binding protein (ADIP) (Afadin DIL domain-interacting protein) (SSX2-interacting protein), SSX2IP, KIAA0923
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SSX2IP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.