



LIMS1 Antibody, FITC conjugated

Product Code	CSB-PA012955EC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48059
Immunogen	Recombinant Human LIM and senescent cell antigen-like-containing domain protein 1 protein (1-325AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM and senescent cell antigen-like-containing domain protein 1 (Particularly interesting new Cys-His protein 1) (PINCH-1) (Renal carcinoma antigen NY-REN-48), LIMS1, PINCH PINCH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LIMS1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.