



CD82 Antibody, FITC conjugated

Product Code	CSB-PA13269C0Rb
Abbreviation	CD82 antigen
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27701
Immunogen	Recombinant Human CD82 antigen protein (111-228AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1) (Suppressor of tumorigenicity 6 protein) (Tetraspanin-27) (Tspan-27) (CD antigen CD82), CD82, KAI1 SAR2 ST6 TSPAN27
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD82