



ERGIC2 Antibody, Biotin conjugated

Product Code	CSB-PA857020ED01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RQ1
Immunogen	Recombinant Human Endoplasmic reticulum-Golgi intermediate compartment protein 2 protein (55-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum-Golgi intermediate compartment protein 2, ERGIC2, ERV41 PTX1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERGIC2