



PRDM10 Antibody, FITC conjugated

Product Code	CSB-PA878858LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NQV6
Immunogen	Recombinant Human PR domain zinc finger protein 10 protein (1-122AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PR domain zinc finger protein 10 (EC 2.1.1.-) (PR domain-containing protein 10) (Tristanin), PRDM10, KIAA1231 PFM7 TRIS
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRDM10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.