



Recombinant Mouse Peroxisome proliferator-activated receptor gamma (Pparg)

Product Code	CSB-EP018424MOb0
Relevance	Nuclear receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the nuclear receptor binds to DNA specific PPAR response elements (PPRE) and modulates the transcription of its target genes, such as acyl-CoA oxidase. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis. ARF6 acts as a key regulator of the tissue-specific adipocyte P2 (aP2) enhancer. Acts as a critical regulator of gut homeostasis by suppressing NF-kappa-B-mediated proinflammatory responses. Plays a role in the regulation of cardiovascular circadian rhythms by regulating the transcription of ARNTL/BMAL1 in the blood vessels
Abbreviation	Recombinant Mouse Pparg protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P37238
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 95% as determined by SDS-PAGE. ≥ 95% as determined by SEC-HPLC.
Sequence	MGETLGDSPVDPEHGAFADALPMSTSQEITMVDTEMPFWPTNFGISSVDLSV MEDHSHSFDIKPFTTVDVFSSISAPHYEDIPFTRADPMVADYKYDLKLQEYQSAI KVEPASPPYYSEKTQLYNRPHEEPSNSLMAIECRVCGDKASGFHYGVHACEG CKGFFRRTIRLKLIDRCDLNCRIHKSRNKCYCRFQKCLAVGMSHNAIRFG RMPQAEKEKLLAEISSDIDQLNPESADLRALAKHLYDSYIKSFPLTKAKARAILT GKTTDKSPFVIYDMNSLMMGEDKIKFKHITPLQEQSKEVAIRIFQGCQFRSVEA VQEITEYAKNIPGFINDLNDQVTLTKYGVHEIYTMLASLMNKDGVLISEGQGF MTREFLKSRLRKPFGDFMEPKFEFAVKFNALELDDSDLAIFIIVILSGDRPGLLN VKPIEDIQDNLLQALELQLKLNHPESQLFAKVLQKMTDLRQIVTEHVQLLHVIK KTETDMSLHPLLQEIYKDLY
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Pparg
Expression Region	1-505aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at



4°C for up to one week.

Tag Info


N-terminal 10xHis-tagged

Mol. Weight


65.1 kDa

Protein Length

Full Length

Image


(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



The purity of Pparg was greater than 95% as determined by SEC-HPLC

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.