




Recombinant uncultured bacterium Flagellin

Product Code	CSB-EP2379CMPe1
Relevance	Flagellin is the subunit protein which polymerizes to form the filaments of bacterial flagella.
Abbreviation	Recombinant Uncultured bacterium Flagellin protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6QA53
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	uncultured bacterium
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVVQHNLQAMNSNRMLGITQKTASKSTEKLSSGYAINRAADNAAGLAISEKMR KQIRGLTQASTNAEDGISSVQTAEGALTEVHMLQRMNELAIQAANGTNSDD RSYIQDEIDQLTQEIDRVAETTKFNETYLLKGDTKNVDAMDYTYSYKAVTTNTV ARASVLAAENTATGMSVSISFAANSGBKVTAADSNNLAKAIRDQGFITTTSTQNG KVYVYGLELNGSDAKANYTVSTVSMEAGTFKILNSNKQVVASVTISTTASFKKVS GMSQIVTAYSVSAAYATGDVYSLYDADGNAISANKLDKYFTAGGATEAGGIAT TLSANSGVPKVYDVLGKEVSAVSIASLTAVKDKTAALKMNFHVGADGTDNN KIKINIEAMTAKSLGVNGLKVGSGSSGTNATNAIEIIAGAIKKVSTQRSALGAVQNR LEHTINNLDNIVENTTAESGIRDTDMATEMVKYSNANILSQAGQSMLAQSNQS NQGVLQLLQ
Research Area	Microbiology
Source	E.coli
Expression Region	1-493aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	51.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.