



ESRP2 Antibody, Biotin conjugated

Product Code	CSB-PA007833LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6T0
Immunogen	Recombinant Human Epithelial splicing regulatory protein 2 protein (615-720AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Epithelial splicing regulatory protein 2 (RNA-binding motif protein 35B) (RNA-binding protein 35B), ESRP2, RBM35B
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	ESRP2