



LETM1 Antibody, FITC conjugated

Product Code	CSB-PA012877LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95202
Immunogen	Recombinant Human Mitochondrial proton/calcium exchanger protein (443-543AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial proton/calcium exchanger protein (Leucine zipper-EF-hand-containing transmembrane protein 1), LETM1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	LETM1