



MZF1 Antibody, HRP conjugated

Product Code	CSB-PA015389LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28698
Immunogen	Recombinant Human Myeloid zinc finger 1 protein (186-354AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myeloid zinc finger 1 (MZF-1) (Zinc finger and SCAN domain-containing protein 6) (Zinc finger protein 42), MZF1, MZF ZNF42 ZSCAN6
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MZF1