



NUP98 Antibody, HRP conjugated

Product Code	CSB-PA016209LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P52948
Immunogen	Recombinant Human Nuclear pore complex protein Nup98-Nup96 protein (965-1110AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear pore complex protein Nup98-Nup96 (EC 3.4.21.-) [Cleaved into: Nuclear pore complex protein Nup98 (98 kDa nucleoporin) (Nucleoporin Nup98) (Nup98); Nuclear pore complex protein Nup96 (96 kDa nucleoporin) (Nucleoporin Nup96) (Nup96)], NUP98, ADAR2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NUP98
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.