



OCLN Antibody, FITC conjugated

Product Code	CSB-PA016263LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q16625
Immunogen	Recombinant Human Occludin protein (285-430AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Occludin, OCLN
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	OCLN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.