



PEAR1 Antibody, HRP conjugated

Product Code	CSB-PA017765LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5VY43
Immunogen	Recombinant Human Platelet endothelial aggregation receptor 1 protein (932-1037AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet endothelial aggregation receptor 1 (hPEAR1) (Multiple epidermal growth factor-like domains protein 12) (Multiple EGF-like domains protein 12), PEAR1, MEGF12
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PEAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.