



PLXNB1 Antibody, HRP conjugated

Product Code	CSB-PA018222LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43157
Immunogen	Recombinant Human Plexin-B1 protein (664-837AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Plexin-B1 (Semaphorin receptor SEP), PLXNB1, KIAA0407 PLXN5 SEP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PLXNB1